

Journal of Chromatography A, 894 (2000) 135-144

JOURNAL OF CHROMATOGRAPHY A

www.elsevier.com/locate/chroma

Development and validation of a capillary electrophoresis method for the measurement of short-chain organic acids in natural rubber latex

V. Galli, N. Olmo, C. Barbas*

Facultad de CC Experimentales y Técnicas, Universidad San Pablo-CEU, Urbanización Montepríncipe, Ctra. Boadilla del Monte km 5,3, 28668 Madrid, Spain

Abstract

Short-chain organic acid contents in serum of natural latex are interesting to measure and capillary electrophoresis (CE) has proved to be a good tool for their study. In the present work a method has been developed to identify the short-chain organic acids present in sera of natural rubber latex (oxalic, formic, fumaric, aconitic, succinic, malic, glutaric, citric, acetic, glycollic, propionic and quinic acids), the separation was optimised and the quantification method validated. The separation was performed on a CE system with UV detection at 200 nm. The separation was carried out with an uncoated fused-silica capillary (57 cm×50 μ m I.D.) and was operated at -10 kV potential. The separation buffers were prepared with 0.5 *M* H₃PO₄, 0.5 m*M* cetyltrimethylammonium bromide and pH adjusted by adding NaOH to 6.25 except for propionic acid which was better measured at pH 7.00. Validation parameters are adequate and limits of detection range from 0.005 m*M* to 1.6 m*M*. Short-chain organic acids were measured with this method in sera of three different types of latex. © 2000 Elsevier Science B.V. All rights reserved.

Keywords: Validation; Rubber; Latex; Organic acids; Fatty acids

1. Introduction

The natural latex concentrates used for the manufacture of dipped goods, such as gloves, are products both of nature and technology. Nature provides the raw material, field latex, from the tree species *Hevea brasiliensis* and it is then modified, preserved, concentrated and tested to produce the commercial concentrates [1].

Nowadays, the increasing use of latex products in healthcare areas has also increased new scientific

controls of raw material and manufactured products in order to maintain higher standards in quality control.

The aqueous phase of fresh natural latex contains many different chemical species, of which the principal types are carbohydrates, electrolytes, proteins and amino acids. In the absence of adequate preservation, the carbohydrates become microbiologically oxidised to the so-called volatile fatty acids (VFAs), which comprise mainly formic, acetic and propionic acids [2]. For this reason their measurement gives an indication of the extent of microorganic activity which has occurred during the time since the latex exuded from the tree. These matters are of consider-

^{*}Corresponding author. Fax: +34-91-3510-475.

E-mail address: cbarbas@ceu.es (C. Barbas).

^{0021-9673/00/\$ –} see front matter © 2000 Elsevier Science B.V. All rights reserved. PII: S0021-9673(00)00525-2

able industrial importance, because several months may elapse between the production of preserved and concentrated natural rubber latex and its eventual use. VFA number measured in the steam of a quantity of latex is a classical test in quality control.

The ionic content of latex is of continuing interest since it is believed to contribute to the properties and especially the processing of the latex, at least for natural rubber latices of low VFA number, it is thought that the salts of non-volatile acids (NVAs) significantly influence latex mechanical stability [3]. According to diverse sources [1,4-6], the principal non-volatile organic acids found in natural rubber latex are citric and malic acids. Minor amounts of trans-aconitic, fumaric, glycollic, succinic, quinic and glutaric acids may also be present. Furthermore, some of these acids mediate the biogenesis of polyisoprenes and their measurement can contribute to the study of the metabolism of the plants which produce them, since their biological function remains obscure.

Previously reported works on the measurement of short-chain organic acids in natural rubber latex were carried out using ion chromatography [1,5]. Capillary electrophoresis (CE) has proved to be a good choice for investigation of samples in aqueous media, since usually no more than a simple dilution of samples is needed. Supporting this statement are the studies of Kenney [7] who reported the separation of some organic acids of interest in a variety of food matrices with indirect detection employing a commercial pack. Huang et al. [8] reported the separation of a six-component mixture of monocarboxylic acids from formate to hexanoate with conductivity detection. Devêvre et al. [9] presented the separation of 14 organic acids with indirect UV detection using a commercial modifier: however, fumaric acid is not detected in the working conditions with this method. Turcat et al. [10] reported the monitoring of six organic acids in snow and rain water and Shirao et al. [11] reported the determination of organic acids in urine. Our research group has developed several methods of separation of organic acids in diverse matrices [12-14] with direct UV detection at 200 nm where the absorbance has a maxim. This low wavelength can be employed when working with non UV-absorbing aqueous buffers.

The aim of the present work was the development of a CE method to identify the short-chain organic acids present in sera of natural rubber latex, optimise the separation and validate the quantification method.

2. Materials and methods

2.1. Instrumentation

The separation was performed on a CE P/ACE 5500 system (Beckman) with UV detection at 200 nm. The injection was by pressure (3.3 bar) for 5 s. The separation was carried out with an uncoated fused-silica capillary (57 cm×50 μ m I.D.) and was operated at -10 kV potential. Temperature was maintained at 25°C. For optimising the pH of the separation buffers were prepared with 0.5 *M* H₃PO₄, 0.5 m*M* cetyltrimethylammonium bromide (CTAB), as cationic surfactant to avoid electroosmotic flow, and pH adjusted by adding NaOH in the range from 5.5 to 7.5. The best pH found for separation was 6.25 except for propionic acid which was better measured at pH 7.00.

2.2. Chemicals

Standards were obtained from Sigma (St. Louis, MO, USA). Phosphoric acid (85%) was from Merck (Darmstadt, Germany), sodium hydroxide from Panreac (Madrid, Spain) and organic solvents from Scharlau (Barcelona, Spain).

2.3. Samples

Samples of Thailand natural latex double centrifuged were directly imported and were representative of commercial material. For valuation, samples of natural latex coming from Guatemala taken directly from the trees or single centrifuged were also used. To obtain serum from latex ultracentrifugation has been utilised, but it needs 30 000 rpm for 60 min, so different chemical procedures for easier colloidal destabilisation were assayed, bearing in mind that many additives could alter or interfere with the electrophoretic process. Therefore, acidic substances such as H_3PO_4 , water-soluble electrolytes as $Ca(NO_3)_2$, $(NH_4)_2SO_4$ or $(NH_4)_3PO_4$ and hydrophobic organic solvents as acetone or chloroform were assayed.

The best conditions found were to mix 10% (v/v)

 H_3PO_4 with latex 1:2 (v/w), to stir vigorously and keep at $-20^{\circ}C$ overnight. The next day, after tempering, serum was filtered through 0.45- μ m nylon filters and diluted 1:3 (v/v) with purified water prior to analysis.

2.4. Validation

Individual stock solutions of each organic acid, 400 mM in water, were prepared and stored at -20° C at this concentration. On the day of the analysis they were adequately diluted.

Linearity of response for standards was tested assaying by triplicate six levels of concentrations, ranging from 0.05 to 0.5 m*M* for nitrate; oxalic, fumaric and aconitic acids; from 0.5 to 10 m*M* for formic, glycollic and quinic acids; from 0.1 to 5.0 m*M* for succinic and glutaric acids; from 0.25 to 2.5 m*M* for malic and citric acids; from 0.5 to 5 m*M* for acetic acid and 0.49 to 7.81 m*M* for propionic acid. Linearity of response for samples was tested in the same way but replacing water with latex serum.

Recovery was estimated comparing the values obtained in the linearity of the calibration serum, with the standards linearity, taking into account the endogenous serum concentrations, which had been previously quantified.

Within-day precision was tested both to check the constancy of instrumental response to a given analyte and the repetitiveness of concentrations and migration times, since the latter is a key parameter for peak assignment. For this purpose, the assay was performed with 10 solutions of standards and 10 of samples, in the medium concentration of the calibration curve for all the compounds.

Limits of detection (LODs) were calculated following IUPAC recommendations $[(a+3S_B)/b]$ for chromatographic methods [15] by extrapolating to zero concentration the standards deviation of the last three points of linearity and interpolating this value in the corresponding equation.

3. Results and discussion

3.1. Coagulation

The addition of acidic substances to latices reduces the colloid stability and is commonly attributed to chemical interaction between the added hydrogen ions and the surface-bound carboxylate ions which confer colloid stability upon the latex [3]. Phosphoric acid was chosen because it is the same acid included in the electrolyte buffer and it is not UV absorbing. The concentration used was varied from 30% to 10% to get the least ionic strength possible in the sample as lower ionic strength in samples than in buffers have been recommended to reduce electrodispersion broadening and extremely high ionic strength in buffers causes high Joule heating.

Different water-soluble electrolytes were also assayed. Monovalent salts need the addition of large amounts to bring about flocculation and this is not useful in the electrophoresis process. Ammonium sulphate and ammonium phosphate were highly effective as destabilising agents, but two big system peaks appeared in the electropherogram overlapping




Fig. 1. Electropherogram of standards at different pH. Peaks: (1) oxalic, (2) formic, (3) fumaric, (4) aconitic, (5) succinic, (6) malic, (7) glutaric, (8) citric, (9) acetic, (10) glycollic, (11) propionic, (12) furanoic and (13) quinic acids. Buffer, 0.5 *M* H₃PO₄-0.5 m*M* CTAB. UV detection at 200 nm. Injection by pressure (3.3 bar) for 5 s. Uncoated fused-silica capillary (57 cm×50 μ m I.D.), -10 kV potential. Temperature 25°C.




Fig. 2. Electropherogram of serum of latex at different pH. For other conditions see the text. (1) nitrate; (2) oxalic, (3) aconitic, (4) succinic, (5) malic, (6) citric, (7) acetic, (8) glycollic and (9) propionic acids; (10) unknown. Buffer, 0.5 M H₃PO₄-0.5 mM CTAB. UV detection at 200 nm. Injection by pressure (3.3 bar) for 5 s. Uncoated fused-silica capillary (57 cm×50 μ m I.D.), -10 kV potential. Temperature 25°C.




Fig. 3. Identification of peaks by spiking samples with standards. Standards: (1) nitrate; (2) oxalic, (3) formic, (4) fumaric, (5) aconitic, (6) succinic, (7) malic, (8) glutaric, (9) citric, (10) acetic, (11) glycollic, (12) propionic and (13) quinic acids. Sample: (1) nitrate; (2) oxalic, (3) formic, (5) aconitic, (6) succinic, (7) malic, (9) citric, (10) acetic and (12) propionic acids; (14) unknown. Buffer, 0.5 M H₃PO₄-0.5 mM CTAB pH 6.25 with NaOH. UV detection at 200 nm. Injection by pressure (3.3 bar) for 5 s. Uncoated fused-silica capillary (57 cm×50 μ m I.D.), -10 kV potential. Temperature 25°C.

with analytes. It is well known that salts of divalent cations such as calcium cause rapid particle agglomeration of the latex and separation of the polymer but, after separation, a white precipitate in aqueous phase was observed, probably due to the calcic salts of the organic acids, and so these salts were avoided.

Finally, none of the assayed organic solvents in the studied range (10-30%, v/w) gave adequate results. Therefore, the conditions established for the method were those explained above with phosphoric acid and freeze-thaw destabilisation to force the process.

3.2. Optimisation

Starting from the experience in our laboratory separating short-chain organic acids in different matrices, described above, we aimed at optimising sensitivity and selectivity of the electrophoretic method by achieving good peak shapes and separation of peaks. For this purpose, the ionic strength was raised by increasing the concentration of the buffer, to get the biggest possible stacking effect. As samples were coagulated in highly concentrated H_3PO_4 the buffer was made 0.5 M in phosphate. It poses a problem with high current intensities and high heat generation by Joule effect, but with the cooling system in the equipment it had no deleterious effects on separation. As some of the short-chain organic acids existing in latex serum partially overlapped under previously developed conditions, pH in the buffer was optimised varying between 5.50 and 7.50 with 0.25 increments. Fig. 1 shows the separation obtained for standards at the different pH values. The resolution for standards is complete in the range 6.00 to 6.50. Under pH 6.00 some peaks disappeared overlapping with others and over pH 6.50 resolution is not complete in some cases. Samples were then analysed in the same range (6.00-6.50) and, as shown in Fig. 2, the best results ocurred at 6.13 and 6.25 so finally the latter value was established. The identification of peaks was performed by comparing with the migration times of pure standards and by spiking the sample with each pure standard separately. Fig. 3 shows the electropherograms under final conditions for standards, a

double centrifuged commercial latex alone and with the working standards added.

The ghost negative peak in the electropherogram is a system peak due to the high electric field created by the elevated phosphoric acid concentration in samples, necessary to coagulate latex. When samples were neutralised with NaOH, previous to injection, at increasing pH, the peak diminished or disappeared but, as could be seen in Fig. 4, measurement is not possible in this area anyway. As this treatment increases the operations without any further improvement, it was avoided.

During identification it could be seen that the system peak appearing in samples made quantification of propionic acid difficult, so a second condition




Fig. 4. Effect on electropherogram profile of neutralisation in samples: (1) nitrate; (2) oxalic, (3) formic, (4) aconitic, (5) succinic, (6) malic, (7) citric, (8) acetic, (9) glycollic and (10) propionic acids; (11) unknown. Buffer, $0.5 M H_3PO_4-0.5 mM$ CTAB. UV detection at 200 nm. Injection by pressure (3.3 bar) for 5 s. Uncoated fused-silica capillary (57 cm×50 µm I.D.), -10 kV potential. Temperature 25°C. A: Sample–water (1:1); B: sample–1% NaOH (1:1); C: sample–2% NaOH (1:1); D: sample–2.5% NaOH (1:1); E: sample–3% NaOH (1:1); F: sample–4% NaOH (1:1); G: sample–5% NaOH (1:1).

was developed just for this acid to displace it out of the negative peak, which was pH 7.00. So for measuring propionic acid a second injection must be done in a different buffer, but everything is done automatically.

An unknown peak appears near 12.5 min. Great efforts have been made in order to identify it. A long list of compounds, described as possible latex constituents and with adequate size for the migration time, has been assayed: aspartic acid, glutaric acid, arginine, alanine, hystidine, glycine, lysine, glycerophosphate, levulinic acid, mevalonic acid, caprylic acid, shikimic acid, isovaleric acid, ketoisovaleric acid, α -ketoglutaric acid, benzoic acid, valproic acid, ftalic acid, lactic acid, maleic acid, methionine, butyric acid, gluthation and ascorbic acid, but currently it remains unknown.

3.3. Validation

As shown in Table 1, standards and samples fit the linear model (r>0.99) for all the organic acids. Although a small bias was found in some of them because the intercept (a)±its limits of confidence (LCs) does not include the zero value, it has statistical significance, but it has no practical consequences as can be seen in the recoveries that are near 100%

Table 1			
Main parameters	of validation	for linearity	and accuracy

		Linearity		Range	Accuracy:	LOD		
		a±LC	b±LC	r	(m <i>M</i>)	Recovery $(\%) \pm RSD(\%)$	(m <i>M</i>)	
Nitrate	Standard	39±79	35 839±889	0.999	0.050-0.500	100±4	0.022	
	Sample	64 ± 94	38 942±1053	0.999	0.016 - 0.166	109±7	0.005	
Oxalic acid	Standard	111 ± 124	9829±615	0.997	0.050 - 0.500	100 ± 12	0.020	
	Sample	-122 ± 141	$12\ 113 \pm 933$	0.992	0.086 - 0.236	104 ± 10	0.071	
Formic acid	Standard	155±216	414±35	0.991	0.500 - 10.00	93±22	1.612	
	Sample	155 ± 139	448±29	0.995	0.375 - 7.500	108±9	1.037	
Fumaric acid	Standard	-81 ± 352	64 898±1735	0.999	0.037-0.375	100±3	0.003	
	Sample	-299 ± 388	67 955±2095	0.999	0.037 - 0.375	100 ± 4	0.003	
Aconitic acid	Standard	-72 ± 504	52 203±2494	0.998	0.050-0.500	99±6	0.002	
	Sample	-81 ± 141	56 856±1595	0.999	0.015 - 0.165	110±5	0.005	
Succinic acid	Standard	-388 ± 332	3977±135	0.998	0.100-5.000	107±12	0.049	
	Sample	-455 ± 172	4165 ± 95	0.999	0.193-3.943	104±9	0.117	
Malic acid	Standard	323±132	4122±131	0.999	0.250-2.500	99±6	0.204	
	Sample	-405 ± 246	5081 ± 235	0.997	0.740 - 1.490	105±6	0.224	
Glutaric acid	Standard	-442 ± 370	4464 ± 150	0.998	0.100 - 5.00	102 ± 8	0.076	
	Sample	-373 ± 338	4808 ± 164	0.999	0.375 - 3.750	113 ± 10	0.046	
Citric acid	Standard	296±193	6234±191	0.999	0.250 - 2.500	100±5	0.159	
	Sample	-205 ± 403	7294±427	0.995	0.635-1.385	107±5	0.140	
Acetic acid	Standard	-286 ± 196	1853±71	0.998	0.500 - 5.000	99±6	0.177	
	Sample	-354 ± 96	1807 ± 51	0.998	0.294 - 4.044	95 ± 10	0.234	
Glycollic acid	Standard	438±437	1184 ± 70	0.995	0.500 - 10.00	99±7	1.414	
	Sample	-187 ± 113	1412 ± 24	0.999	0.375 - 7.500	102 ± 8	0.185	
Propionic acid	Standard	-135 ± 173	889±36	0.998	0.488-7.805	100 ± 2	0.650	
	Sample	320±137	844±38	0.998	0.366-5.854	111 ± 8	1.600	
Quinic acid	Standard	1011 ± 1138	3007±182	0.995	0.500-10.00	96±12	1.188	
	Sample	113 ± 213	3593 ± 45	0.999	0.375 - 7.500	101 ± 15	0.413	

in the whole range and it can be due to the good fit of the points to the regression line. During the calculation of recoveries it could be detected that for fumaric acid a recovery of 140% was obtained if pH of standards was not the same as in samples, so standards were prepared with the same quantity of 10% (v/v) H_3PO_4 as in samples and all the recoveries obtained were near 100%.

Mathematical LODs range from 0.005 mM to 1.6 mM. It is a problem, due to the statistical approach, that some of them are over the validated range.

Obviously, the method is applicable for quantification in the range experimentally proved, which many times is under the calculated detection limit. It is interesting to see the doubtful meaning of a mathematical approach. In the same line, samples have lower detection limits than calculated for standards in some acids. That is justified because when the sample contents one of the acids in high concentration the lowest points in the regression line are not so low, and the corresponding standard deviation are smaller, giving better detection limits.

Table 2			
Main validation	parameters	for	precision

		Precision							
		Intra-assay			Inter-assay				
		Mean (m <i>M</i>)	RSD (%)	t _m ^a (min)	RSD (%)	Mean (m <i>M</i>)	RSD (%)	t _m ^a (min)	RSD (%)
Nitrate	Standard	0.244	2	7.162	0.04	0.244	6	7.080	1
	Sample	0.178	2	7.192	0.20	0.188	8	7.094	1
Oxalic acid	Standard	0.244	5	7.824	0.45	0.244	9	7.874	1
	Sample	0.341	4	7.857	0.38	0.319	9	7.914	1
Formic acid	Standard	4.879	5	8.213	0.37	4.879	7	8.108	2
	Sample	2.771	6	8.262	0.21	2.820	8	8.130	1
Fumaric acid	Standard	0.244	2	8.772	0.31	0.244	6	8.634	2
	Sample	0.131	1	8.830	0.23	0.136	6	8.681	1
Aconitic acid	Standard	0.244	2	8.927	0.83	0.244	8	8.814	2
	Sample	0.162	1	9.027	0.21	0.170	7	8.869	1
Succinic acid	Standard	2.444	1	9.247	0.37	2.444	4	9.115	2
puttine und	Sample	1.510	2	9.352	0.24	1.572	5	9.281	1
Malic acid	Standard	1.220	3	9.442	0.37	1.220	4	9.304	2
interio dere	Sample	2.276	2	9.527	0.26	2.319	4	9.355	1
Glutaric acid	Standard	2.444	2	9.892	0.83	2.444	6	9.582	3
	Sample	1.264	2	9.727	0.26	1.306	7	9.547	1
Citric acid	Standard	1.220	1	9.918	0.38	1.220	5	9.771	2
chine uelu	Sample	2.383	2	9.990	0.27	2.473	5	9.799	1
Acetic acid	Standard	2.444	4	10.110	0.38	2.444	7	9.953	2
ricelle della	Sample	1.620	4	10.175	0.26	1.664	5	9.978	2
Glycollic acid	Standard	4.878	2	10.550	0.39	4.878	4	10.382	2
Sigeonic acid	Sample	2.251	2	10.558	0.25	2.418	6	10.345	2
Propionic acid	Standard	3.902	2	10.553	0.36	3.902	3	10.485	1
r ropionie aciu	Sample	2.394	6	10.555	0.30	2.591	8	10.485	2
Quinic acid	Standard	4.878	3	13.963	0.51	4.878	4	13.677	3
Quille aciu	Sample	2.598	3	13.903	0.31	2.772	4 5	13.780	2

^a t_m = Migration time.

	Anion concentration (mmol l^{-1} serum)					
	A	В	С	D	Е	
Nitrate	0.14	n.d	n.d	_	_	
Oxalic acid	0.82	0.52	0.37	0.97 - 1.99	0.12-1.37	
Formic acid	3.13	9.50	6.92	2.00-6.66	n.d-8.22	
Fumaric acid	n.d	n.d	n.d	_	_	
Aconitic acid	0.14	0.23	0.16	_	_	
Succinic acid	1.51	23.10	16.68	n.d-5.00	n.d-10.25	
Malic acid	6.27	1.80	1.90	6.39-11.66	4.21-10.14	
Glutaric acid	n.d	n.d	n.d	_	_	
Citric acid	7.36	8.08	5.74	2.00-7.73	3.73-6.89	
Acetic acid	2.95	24.50	17.04	6.39-16.77	n.d-18.15	
Glycollic acid	n.d	n.d	n.d	_	-	
Propionic acid	2.83	n.d	n.d	_	_	
Quinic acid	n.d	n.d	n.d	_	_	

Table 3 Short-chain organic acid content in sera of different types of latices^a

^a Sample A: latex double centrifuged from Thailand (DRC 60.27%); sample B: natural latex from Guatenmala (DRC 29/85%); sample C: natural latex single centrifuged from Guatemala (DRC 58.46%); sample D: high ammoniated concentrated latex from Malaysia (DRC 59.73%) [5]; sample E: low ammoniated concentrated latex from Malaysia (DRC 59.78%) [5].

When running 10 runs per day of both standards and samples, daily relative standard deviation (RSDs) in concentrations are low enough to consider the method acceptable (1 to 6%) and variations when they exist are mainly due to difficulties in integration of small peaks. RSDs for migration times, Table 2, are very low when they are measured in the same assay for both standards and samples (0.04 to 0.8%).

The intermediate precision evaluated on different days with a total of 25 runs provided RSD values




Fig. 5. Electropherogram of serum obtained by freezing. Peaks: (1) nitrate; (2) thiosulfate; (3) oxalic, (4) formic, (5) aconitic, (6) succinic, (7) malic, (8) citric and (9) acetic acids; (10) unknown. Buffer, 0.5 M H₃PO₄-0.5 mM CTAB. UV detection at 200 nm. Injection by pressure (3.3 bar) for 5 s. Uncoated fused-silica capillary (57 cm×50 μ m I.D.), -10 kV potential. Temperature 25°C.

slightly superior to intra-assay precision (3 to 9% for concentrations and 1 to 3% for migration times), as could be expected.

3.4. Quantification of samples

After development, the method was applied to the quantification of these acids in serum of three different types of latices: the one double centrifuged coming from Thailand used during validation (sample A), one just obtained from the tree without any treatment (sample B) and one centrifuged once (sample C), both coming from Guatemala. Results are shown in Table 3. As it could be expected, natural latex without any treatment has higher organic acids content in the same weigh, mainly because during the centrifugation process the ratio of rubber increases while the serum is partially eliminated. In these acids DRC (dry rubber content) is 60.27% for A. 29.85% for B and 58.46% for C. But. variations in this parameter do not fully explain differences in short-chain organic acids content. So further studies must be done to establish the correlation between these acids and different processes in latex. When comparing the values obtained and those existing in the literature [1,5,6] it can be found that the published ranges are wide and so we are in agreement with them in most cases, but it is not very clarifying.

3.5. Coagulation by freezing

While assaying the best way to receive samples from different countries to apply the method in phytochemical studies, the possibility of destabilising the latex by freezing at -20° C overnight, without any other addition could be seen. This is very interesting because, just the acids naturally free in serum are measured, avoiding the liberation of acids bound to solid matter in some way or those lost due to insolubilisation of some acids in high acidic media. The method consisted of filling 50-ml Falcon plastic tubes with around 37 g of latex and placing them in the freezer at -20° C overnight, after tempering, the tubes were centrifuged at 4000 g for 20 min and serum was filtered and injected. As the ionic strength is lower than in the previous method, favouring stacking effect, and no dilution is needed, Table 4

Results obtained from 10 samples coagulated with phosphoric acid	l
and 10 samples coagulated by freezing ^a	

	10% H_3PO_4 (mmol 1^{-1} serum)	Frozen (mmol 1^{-1} serum)
Nitrate	0.18 ± 0.02	0.76 ± 0.02
Oxalic acid	0.91 ± 0.06	1.05 ± 0.02
Aconitic acid	$0.14 {\pm} 0.01$	0.16 ± 0.01
Succinic acid	1.24 ± 0.12	1.41 ± 0.03
Malic acid	8.00 ± 0.41	8.85 ± 0.20
Citric acid	6.25 ± 0.29	7.01 ± 0.23
Acetic acid	2.65 ± 0.31	3.04±0.19

^a Mean \pm limits of confidence (P > 95%).

the electropherogram profile improved considerably, as can be observed in Fig. 5.

Only a new peak appeared in the electropherogram due to an unfolding of the nitrate peak previously missed. It was identified as tiosulfate by retention time and by spiking the sample with it.

In order to compare the results obtained coagulating with phosphoric acid and just freezing, two series of 10 samples of the same latex were processed individually coagulating by both methods. Results can be found in Table 4. Only nitrate poses a great difference, which was previously justified. Concentration of organic acids are slightly superior and limits of confidence slightly inferior in the frozen way than when coagulating with phosphoric acid. The value could be justified because losses by volatisation are minimum and the smaller deviations because the peaks are better measured.

4. Conclusions

A CE method has been developed and validated for measuring short-chain organic acids in sera of natural rubber latex directly without any sample pre-treatment more than coagulation. It has been applied to different types of samples but further studies must be done to establish the correlation between these acids and different processes in latex.

Acknowledgements

The present work has been developed with the collaboration of Tecnilatex S.A.

References

- [1] T.D. Pendle, Rubber Chem. Technol. 63 (3) (1990) 234.
- [2] D.C. Blackley, Polymer Latices Science and Technology, 2nd ed., Types of Latices, Vol. 2, Chapman and Hall, London, 1997.
- [3] D.C. Blackley, Polymer Latices Science and Technology, 2nd ed., Fundamental Principles, Vol. 1, Chapman and Hall, London, 1997.
- [4] K.C. Calvert, Plastics Rubber: Mater. Appl. 2 (1977) 59.
- [5] R.C. Crafts, A.D.T. Gorton, T.D. Pendle, NR Technol. 16 (1985) 12.
- [6] R.T. Davies, T.D. Pendle, Rubber Dev. 44 (4) (1991) 94.
- [7] B.F. Kenney, J. Chromatogr. 546 (1991) 423.

- [8] X. Huang, J.A. Luckey, M.J. Gordon, R.N. Zare, Anal. Chem. 61 (1989) 766.
- [9] O. Devèvre, D.P. Putra, B. Botton, J. Garbaye, J. Chromatogr. A 679 (1994) 349.
- [10] S. Turcat, P. Masclet, T. Lissol, Sci. Total Eviron. 158 (1994) 21.
- [11] M. Shirao, R. Furuta, S. Suzuki, H. Nakazawa, S. Fujita, T. Maruyama, J. Chromatogr. A 680 (1994) 247.
- [12] C. Barbas, N. Adeva, R. Aguilar, M. Rosillo, T. Rubio, M. Castro, Clin. Chem. 44 (1998) 1340.
- [13] A. García, C. Barbas, R. Aguilar, M. Castro, Clin. Chem. 44 (1998) 1905.
- [14] C. Barbas, J.A. Lucas, F.J. Gutiérrez-Mañero, Phytochem. Anal. 10 (1999) 55.
- [15] G.L. Long, J.D. Winefordner, Anal. Chem. 55 (1983) 712.